



## **Thought Technology Ltd.**

2180 Belgrave Avenue, Montreal Quebec, H4A 2L8 Canada  
Tel: 800-361-3651 (514) 489-8251 Fax: (514) 489-8255  
Email: [mail@thoughttechnology.com](mailto:mail@thoughttechnology.com)  
[www.thoughttechnology.com](http://www.thoughttechnology.com)

### **BioGraph Infiniti 5.1 Computer Requirements:**

#### **Recommended**

- Desktop or laptop with two monitor capacity.
- CPU: 1.8GHz dual core processor or faster.
- Operating systems:
  - Microsoft® Windows® 7 Home Premium, Business, Ultimate, or Enterprise
  - Microsoft® Windows Vista® Home Basic, Home Premium, Business, Ultimate, or Enterprise with Service Pack 2
  - Microsoft® Windows® XP Professional or Home Edition with Service Pack 3
- DirectX 9c or later (if this is missing, it will be installed with BioGraph Infiniti).
- 50 – 60 gigabytes hard disk space for video recording and processing. (The software needs 2.5 gigabytes to install and run.)
- Minimum RAM: 1 GB (Windows® XP) or 2 GB (Windows Vista® or Windows® 7).
- CD ROM or DVD drive.
- Video Card (minimum XGA, 1,024 x 768 monitor resolution)
- 32 bit compatible sound card & speakers
- 1 to 4 USB ports, depending on the desired number of Pro/Flex Infiniti encoders
- 1 to 4 serial com ports, depending on the desired number of ProComp+ encoders.
- Mouse or compatible pointing device
- Microsoft® Office with Word® Excel® (for report generation and printing)
- Compact Flash Reader (for use with compact flash card only)
- Webcam 30 frames per second (for video purposes only)
- Internet access (for updating software and receiving technical support)

#### **PLEASE NOTE:**

- **Tele-Infiniti CF version 3.1 (T9600) is not compatible with 64-bit Windows® operating systems.**
- *When using Windows® 7 Home Basic or Business or Windows Vista® Home Basic or Business, additional codecs must be installed to view camera and DVD.*
- *Low speed CPUs (1.8 GB or less) may be sluggish when running screens with multiple instruments.*